

## Khione Dipper Ice-Bath Setup Manual

Congratulations on your purchase of the Khione Dipper Ice-Bath! With its portable design and quick setup, you'll be enjoying the benefits of cold therapy in no time. Here's a simple guide to setting up your Khione Dipper.

### What's in the Box?

- 1 Khione Dipper Ice-Bath Tub
- 8 Support Poles
- 1 Top Cover
- 1 Air Pump
- 1 Drainage Valve with Rubber Seal

### Setup Instructions:

#### 1. Prepare for Assembly


- Find a flat, clean surface where you plan to use the Khione Dipper.
- Ensure the area is free of sharp objects to avoid damaging the tub.

#### 2. Unpack the Contents

- Carefully remove all parts from the package.
- Lay out the 8 poles, top ring, top cover, and the air pump for easy access.

#### 3. Assemble the Frame

- Locate the 8 entrances at the base of the tub.
- Insert each of the 8 poles into these entrances to form the support frame of the Dipper.



#### 4. Inflate the Top Ring

- Attach the provided air pump to the valve on the top ring.
- Pump to inflate the ring until it is fully expanded and firm.

#### 5. Inflate the Top Cover

- Similarly, attach the pump to the valve on the top cover.
- Inflate until fully expanded.

## 6. Attach the Drainage Valve

- Turn the tub upside down to access the bottom.
- Place the drainage valve onto the designated area, making sure the text on the valve is facing outward.
- Secure the valve by using the provided rubber ring to ensure a tight seal.



## 7. Position Your Dipper

- Carefully turn the tub right side up and place it in the desired location.

## 8. Fill with Water

- Fill the Khione Dipper with water to the recommended level.
- Adjust the water temperature to your preference for cold therapy.

## 9. Enjoy

- Once filled, your Khione Dipper is ready for use.
- Enter carefully and enjoy the refreshing benefits of your very own ice-bath!

### Important Tips:

- Always check the drainage valve for a secure fit before each use to prevent leaks.
- When not in use, store the Khione Dipper in a cool, dry place away from direct sunlight.

Thank you for choosing the Khione Dipper Ice-Bath. For any questions or further assistance, please visit our website at [www.khione.store](http://www.khione.store) or contact our customer support team.

Enjoy your cold therapy journey with the Khione Dipper!



